

# LAMPIRAN

## 1. Code App.py

```
import streamlit as st
from detect import *
from PIL import Image #membuka atau akses gambar
import cv2
import time
import psutil
import numpy as np
import argparse
import os
import platform
import sys
from pathlib import Path
from collections import Counter
import torch.backends.cudnn as cudnn
import torch
from numpy import random
from utils.google_utils import gsutil_getsize
from utils.metrics import fitness
from utils.torch_utils import init_torch_seeds
from models.experimental import attempt_load
from utils.datasets import LoadStreams, LoadImages
from utils.general import check_img_size, check_requirements, check_imshow, non_max_suppression,
apply_classifier, \
    scale_coords, xyxy2xywh, strip_optimizer, set_logging, increment_path
from utils.plots import plot_one_box
from utils.torch_utils import select_device, load_classifier, time_synchronized, TracedModel

#-----Web Page Designing-----
hide_menu_style = """
<style>
    MainMenu {visibility: hidden;}
    div[data-testid="stHorizontalBlock"]> div:nth-child(1)
    {
        border : 2px solid #d0e0db;
        border-radius:5px;
        text-align:center;
        color:black;
        background:dodgerblue;
        font-weight:bold;
        padding: 25px;
    }
    div[data-testid="stHorizontalBlock"]> div:nth-child(2)
    {
        border : 2px solid #d0e0db;
        background:dodgerblue;
        border-radius:5px;
        text-align:center;
        font-weight:bold;
        color:black;
        padding: 25px;
    }
</style>
"""
main_title = """
<div>
    <h1 style="color:dodgerblue;
    text-align:center; font-size:35px;
    margin-top:-95px;">
    CALCULATED POLLEN VIABILITY</h1>
</div>
"""
sub_title = """
<div>
    <h6 style="color:dodgerblue;
    text-align:center;
    margin-top:-40px;">
    Polen Object Detection and Counting </h6>
</div>
"""
#-----Main Function for Execution-----
def main():
    st.set_page_config(page_title='Dashboard',
        layout = 'wide',
        initial_sidebar_state = 'auto')
    st.markdown(hide_menu_style,
        unsafe_allow_html=True)
```

```

st.markdown(main_title,
             unsafe_allow_html=True)
st.markdown(sub_title,
             unsafe_allow_html=True)
inference_msg = st.empty()
st.sidebar.title("Heart Made wkwkwkwkwk")
activities = ["Polen", "Kecambah"] #DIUBAH TANGGAL 23/02/2023
choice = st.sidebar.selectbox("Please Select Activity", activities)
if choice == "Polen":
    input_source = 'Image'
    conf_thres = st.sidebar.text_input("Class confidence threshold",
                                       "0.5")

    save_output_image = st.sidebar.radio("Save output?",
                                       ('Yes', 'No'))


    if save_output_image == 'Yes':
        nosave = False
        display_labels=False
    else:
        nosave = True
        display_labels = True
    weights = "yolov7 (1).pt"
# ----- LOCAL IMAGE -----
if input_source == "Image":

    image = st.sidebar.file_uploader("Select input",
                                     type=["jpg", "jpeg", "png"],
                                     accept_multiple_files=False)

    # nama_file=image.name
    if image is not None:
        our_image = Image.open(image)
        our_image.save("data/gambaranya/{}".format(image.name))
        st.image(our_image)
    if st.sidebar.button("Start Counting"):
        stframe = st.empty()
        st.markdown("""<h4 style="color:black;">
                    Memory Overall Statistics</h4>""",
                    unsafe_allow_html=True)
        kpi5, kpi6 = st.columns(2)
        with kpi5:
            st.markdown("""<h5 style="color:black;">
                        CPU Utilization</h5>""",
                        unsafe_allow_html=True)
            kpi5_text = st.markdown("0")
        with kpi6:
            st.markdown("""<h5 style="color:black;">
                        Memory Usage</h5>""",
                        unsafe_allow_html=True)
            kpi6_text = st.markdown("0")
        a = detect(weights='best.pt',
                  source="data/gambaranya/{}".format(image.name),
                  stframe=stframe,
                  #kpi5_text=kpi5_text,
                  #kpi6_text = kpi6_text,
                  conf_thres=float(conf_thres),
                  device="cpu",
                  classes=[0,1],nosave=nosave,
                  name='{}'.format(image.name)
                  )
        inference_msg.success("Inference Complete!")
        st.image(a, channels="BGR",use_column_width=True)
        torch.cuda.empty_cache()
# ----- MAIN FUNCTION CODE -----
# print(f'Done.({time.time() - t0:.3f}s)')
if __name__ == "__main__":
    try:
        main()
    except SystemExit:
        pass

```

## 2. Code Detect.py



```
1 import argparse
2 import time
3 from pathlib import Path
4
5 import cv2
6 import torch
7 import torch.backends.cudnn as cudnn
8 from numpy import random
9
10 from models.experimental import attempt_load
11 from utils.datasets import LoadStreams, LoadImages
12 from utils.general import check_img_size, check_requirements, check_imshow, \
13     non_max_suppression, apply_classifier, \
14     scale_coords, xyxy2xywh, strip_optimizer, set_logging, increment_path
15 from utils.plots import plot_one_box
16 from utils.torch_utils import select_device, load_classifier, time_synchronized, \
17     TracedModel
18
19 def detect(source='data/gambaranya', img_size=640, conf_thres=0.25, iou_thres=0.45,
20     device='cpu', view_img=False, save_txt=False, save_conv=False, nosave=False,
21     classes="",
22     agnostic_nms=True, augment=False, update=False, project='detect',
23     name='exp', exist_ok=False, no_trace=False, weights='best.pt',
24     stframe=None, save_img=False):
25     source, weights, view_img, save_txt, imgsz, trace = str(source), 'best.pt',
26         True, False, 640, True
27     save_img = True and not source.endswith('.txt') # save inference images
28     webcam = source.isnumeric() or source.endswith('.txt') or
29         source.lower().startswith(
30             ('rtsp://', 'rtmp://', 'http://', 'https://'))
31
32     # Directories
33     save_dir = Path(increment_path(Path('detect') / 'exp', exist_ok=False)) #
34         increment run
35     (save_dir / 'labels' if save_txt else save_dir).mkdir(parents=True,
36         exist_ok=True) # make dir
37
38     # Initialize
39     set_logging()
40     device = select_device('cpu')
41     half = device.type != 'cpu' # half precision only supported on CUDA
42
43     # Load model
44     model = attempt_load(weights, map_location=device) # load FP32 model
45     stride = int(model.stride.max()) # model stride
46     imgsz = check_img_size(imgsz, s=stride) # check img_size
47
48     if trace:
49         model = TracedModel(model, device, 640) #in
50
51     if half:
52         model.half() # to FP16
```

```

45 # Second-stage classifier
46 classify = False
47 if classify:
48     modelc = load_classifier(name='resnet101', n=2) # initialize
49     modelc.load_state_dict(torch.load('weights/resnet101.pt',
50                                     map_location=device)['model']).to(device).eval()
51
52 # Set Dataloader
53 vid_path, vid_writer = None, None
54 if webcam:
55     view_img = check_imshow()
56     cudnn.benchmark = True # set True to speed up constant image size inference
57     dataset = LoadStreams(source, img_size=imgsz, stride=stride)
58 else:
59     dataset = LoadImages(source, img_size=imgsz, stride=stride)
60
61 # Get names and colors
62 names = model.module.names if hasattr(model, 'module') else model.names
63 # colors = [[random.randint(0, 255) for _ in range(3)] for _ in names]
64 colors = [[0,255,0], [0,0,255]]
65
66 # print(colors)
67
68 # Run inference
69 if device.type != 'cpu':
70     model(torch.zeros(1, 3, imgsz,
71                       imgsz).to(device).type_as(next(model.parameters()))) # run once
72 old_img_w = old_img_h = imgsz
73 old_img_b = 1
74
75 t0 = time.time()
76 for path, img, im0s, vid_cap in dataset:
77     img = torch.from_numpy(img).to(device)
78     img = img.half() if half else img.float() # uint8 to fp16/32
79     img /= 255.0 # 0 - 255 to 0.0 - 1.0
80     if img.ndimension() == 3:
81         img = img.unsqueeze(0)
82
83     # Warmup
84     if device.type != 'cpu' and (old_img_b != img.shape[0] or old_img_h !=
85                                img.shape[2] or old_img_w != img.shape[3]):
86         old_img_b = img.shape[0]
87         old_img_h = img.shape[2]
88         old_img_w = img.shape[3]
89         for i in range(3):
90             model(img, augment=False)[0]
91
92     # Inference
93     t1 = time_synchronized()
94     with torch.no_grad(): # Calculating gradients would cause a GPU memory
95         leak
96         pred = model(img, augment=False)[0]

```

```

93     t2 = time_synchronized()
94
95     # Apply NMS
96     pred = non_max_suppression(pred, 0.25, 0.45, classes=[0, 1], agnostic=False)
97     t3 = time_synchronized()
98
99     # Apply Classifier
100    if classify:
101        pred = apply_classifier(pred, modelc, img, im0s)
102
103    # Process detections
104    for i, det in enumerate(pred): # detections per image
105        if webcam: # batch_size >= 1
106            p, s, im0, frame = path[i], '%g: ' % i, im0s[i].copy(),
107                dataset.count
108        else:
109            p, s, im0, frame = path, '', im0s, getattr(dataset, 'frame', 0)
110
111        p = Path(p) # to Path
112        save_path = str(save_dir / p.name) # img.jpg
113        txt_path = str(save_dir / 'labels' / p.stem) + (' ' if dataset.mode ==
114            'image' else f'_{frame}') # img.txt
115        gn = torch.tensor(im0.shape)[1, 0, 1, 0] # normalization gain whwh
116        if len(det):
117            # Rescale boxes from img_size to im0 size
118            det[:, :4] = scale_coords(img.shape[2:], det[:, :4],
119                im0.shape).round()
120
121            # Print results
122            for c in det[:, -1].unique():
123                n = (det[:, -1] == c).sum() # detections per class
124                s += f"{n} {names[int(c)]}'s' * (n > 1)}, " # add to string
125
126            # Write results
127            for *xyxy, conf, cls in reversed(det):
128                if save_txt: # Write to file
129                    xywh = (xyxy2xywh(torch.tensor(xyxy).view(1, 4)) /
130                        gn).view(-1).tolist() # normalized xywh
131                    line = (cls, *xywh, conf) if False else (cls, *xywh) #
132                        label format
133                    with open(txt_path + '.txt', 'a') as f:
134                        f.write((' %g ' * len(line)).rstrip() % line + '\n')
135
136                if save_img or view_img: # Add bbox to image
137                    label = f'{names[int(cls)]} {conf:.2f}'
138                    plot_one_box(xyxy, im0, label=label,
139                        color=colors[int(cls)], line_thickness=1)

```

```

139     j = s.split(" ")
140     # print(j)
141     hidup = []
142     mati = []
143
144     if j[2] != "":
145         mati.append(j[2])
146         hidup.append(j[0])
147         im0 = cv2.putText(im0, str("jumlahnya : "+ str(hidup[0]) + " yang
             hidup, "+ str(mati[0])+" yang mati,"+ " "), (100,100),
             cv2.FONT_HERSHEY_SIMPLEX, 1, (0,0,0), 2, cv2.LINE_AA)
148         # im0 = cv2.putText(im0, str("viabilitasnya : "+ str(hidup) + " yang
             hidup, "+ str(mati)+" yang mati,"+ " "), (100,100),
             cv2.FONT_HERSHEY_SIMPLEX, 1, (0,0,0), 2, cv2.LINE_AA)
149     else:
150         hidup.append(j[0])
151         im0 = cv2.putText(im0, str("jumlahnya : "+ str(hidup[0]) + " yang
             hidup, "), (100,100), cv2.FONT_HERSHEY_SIMPLEX, 1, (0,0,0), 2,
             cv2.LINE_AA)
152     # Stream results
153     if view_img:
154         cv2.imshow(str(p), im0)
155         cv2.waitKey(1) # 1 millisecond
156
157     # Save results (image with detections)
158     if save_img:
159         if dataset.mode == 'image':
160             cv2.imwrite(save_path, im0)
161             print(f" The image with the result is saved in: {save_path}")
162         else: # 'video' or 'stream'
163             if vid_path != save_path: # new video
164                 vid_path = save_path
165             if isinstance(vid_writer, cv2.VideoWriter):
166                 vid_writer.release() # release previous video writer
167             if vid_cap: # video
168                 fps = vid_cap.get(cv2.CAP_PROP_FPS)
169                 w = int(vid_cap.get(cv2.CAP_PROP_FRAME_WIDTH))
170                 h = int(vid_cap.get(cv2.CAP_PROP_FRAME_HEIGHT))
171             else: # stream
172                 fps, w, h = 30, im0.shape[1], im0.shape[0]
173                 save_path += '.mp4'
174             vid_writer = cv2.VideoWriter(save_path,
                 cv2.VideoWriter_fourcc(*'mp4v'), fps, (w, h))
175             vid_writer.write(im0)
176
177     if save_txt or save_img:
178         s = f"\n{len(list(save_dir.glob('labels/*.txt')))} labels saved to
             {save_dir / 'labels'}" if save_txt else ''
179         #print(f"Results saved to {save_dir}{s}")
180
181     print(f'Done. ({time.time() - t0:.3f}s){s}')
182     return im0

```